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## 10.3 FCC SDoC 证书

### Verification of Compliance

This Verification of Compliance is granted to below applicant for following designated equipment

Verification No.:	SZEM2108007269ATV
Applicant:	Shenzhen Yuejiang Technology Co., Ltd.
Address of Applicant:	Room 1003, Building 2, Chongwen Park, Nanshan iPark, No.3370, Liuxian Blvd, Fuguang Community, Taoyuan Street, Nanshan District, Shenzhen
Manufacturer:	Shenzhen Yuejiang Technology Co., Ltd.
Address of Manufacturer:	Room 1003, Building 2, Chongwen Park, Nanshan iPark, No.3370, Liuxian Blvd, Fuguang Community, Taoyuan Street, Nanshan District, Shenzhen
Factory:	Dongguan Xinyou Intelligent Technology Co., Ltd.
Address of Factory:	No.1 Xinyang Road, Dengwu village, Qiaotoutown, Dongguan City
Product Name:	DOBOT MG400/DOBOT Magician Pro
Model No.:	DT-MG400-4R075-01, DT-MGPro-4R075-01
Trade Mark:	DOBOT
Test Report Number(s):	SZEM210800726901

based on tests conducted by SGS-CSTC Standards Technical Services Co., Ltd. Shenzhen Branch on submitted samples of above mentioned product and found to comply with the technical requirements set out in ANSI C63.4 & 47 CFR PART 15B regulations for the evaluation of electromagnetic compatibility.

**Please note:**

- (1) Applicant and/or manufacturer shall fully fulfill FCC SDoC authorization procedure set out in 47 CFR §2.906, §2.909, §2.935, §2.1074, §2.1077.
- (2) Applicant and/or manufacturer shall comply additional requirements set out in §15.21, §15.105 and §15.19 to fully comply regulation 47 CFR Part 15.
- (3) This VoC is ONLY a conclusion of test result, ONLY serves as technical specifications compliance verification for aforementioned product.
- (4) This verification is only valid in conjunction with test report(s) detailed above.

Keny Xu  
EMC Laboratory Manager  
Issue Date:2021-09-02

Copyright of this verification is owned by SGS-CSTC Standards Technical Services Co., Ltd. and may not be reproduced other than in full and with the prior approval of the General Manager. This verification is subjected to the governance of the General Conditions of Services which can accessible at <https://www.sgs.com/en/terms-and-conditions>

Member of SGS Group (Société Générale de Surveillance)

EMC-VOC-F01/ Rev.1.0/ 2020-11-05

## 10.4 CE-MD 证书

**VERIFICATION OF MD COMPLIANCE**

No.:	SZES2101000040MDC
Applicant:	Shenzhen Yuejiang Technology Co., Ltd. Room 1003, Building 2, Chongwen Park, Nanshan iPark, No.3370, Liuxian Boulevard, Fuguang Community, Taoyuan Street, Nanshan District, Shenzhen, Guangdong, China
Manufacturer:	Same as applicant
Product Description:	Industrial Robot
Model No.:	DT-MG400-4R075-01
Sufficient samples of the product have been tested and found to be in conformity with	
Test Standard:	EN ISO 12100:2010 IEC 60204-1:2016 / EN 60204-1:2018 EN ISO 10218-1:2011
as shown in the Test Report Number(s)	SZES210100004001 SZES210110004001

This Verification of MD Compliance has been granted to the applicant based on the results of tests performed by Laboratory of SGS-CSTC Standards Technical Services Co., Ltd. Shenzhen Branch on sample of the above-mentioned product in accordance with the provisions of the relevant specific standards and the Machinery Directive 2006/42/EC. The CE mark as shown below can be affixed, under the responsibility of the manufacturer, after completion of an EC Declaration of Conformity and compliance with all relevant EC Directives. The affixing of the CE marking presumes in addition that the conditions in annexes III of the Directive are fulfilled.

  
 Rocky Wang  
 Laboratory Manager  
 SGS-CSTC

2021-03-12

Copyright of this verification is owned by SGS-CSTC Standards Technical Services (Shenzhen) Co., Ltd. and may not be reproduced other than in full and with the prior approval of the General Manager. This verification is subjected to the governance of the General Conditions of Services, printed overleaf.

Member of SGS Group (Société Générale de Surveillance)

SZE-VOC-F013 / V1.0  
Effective Date : 2020-02-26

This document is issued, on the Client's behalf, by the Company under its General Conditions of Service printed overleaf. The Client's attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein.

Any other holder of this document is advised that information contained herein reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not exonerate parties to a transaction from exercising all their rights and obligations under the transaction documents.

SGSPAPER  
 19947958

## 10.5 CE-EMC 证书

### VERIFICATION OF EMC COMPLIANCE

Verification No.: SZEM2101000072ATV  
Applicant: Shenzhen Yuejiang Technology Co., Ltd  
Address of Applicant: Room 1003, Building 2, Chongwen Park, Nanshan iPark, No.3370, Liuxian Blvd, Fuguang Community, Taoyuan Street, Nanshan District, Shenzhen  
Manufacturer: Shenzhen Yuejiang Technology Co., Ltd  
Address of Manufacturer: Room 1003, Building 2, Chongwen Park, Nanshan iPark, No.3370, Liuxian Blvd, Fuguang Community, Taoyuan Street, Nanshan District, Shenzhen  
Factory: Dongguan Xinyou Intelligent Technology Co., Ltd  
Address of Factory: No.1 Xinyang Road, dengwu village, Qiaotou town, Dongguan City  
Product Description: DOBOT MG400  
Model No.: DT-MG400-4R075-01  
Sufficient samples of the product have been tested and found to be in conformity with  
Test Standards: EN IEC 61000-6-2: 2019  
EN IEC 61000-6-4: 2019  
As shown in the  
Test Report Number(s): SZEM210100007201

This verification of EMC Compliance has been granted to the applicant based on the results of the tests, performed by laboratory of SGS-CSTC Standards Technical Services Co., Ltd. on the sample of the above-mentioned product in accordance with the provisions of the relevant specific standards under Directive 2014/30/EU.

The CE mark as shown below can be used, under the responsibility of the manufacturer, after completion of an EU Declaration of Conformity and compliance with all relevant EU Directives.

Kerry Xu  
EMC Laboratory Manager

Date: 2021-03-12

Copyright of this verification is owned by SGS-CSTC Standards Technical Services Co., Ltd. and may not be reproduced other than in full and with the prior approval of the General Manager. This verification is subjected to the governance of the General Conditions of Services which can be accessible at: <https://www.sgs.com/en/terms-and-conditions>

Member of SGS Group (Société Générale de Surveillance)

EMC-VOC-F01/ Rev.1.0/ 2020-11-05

## 10.6 RoHS 证书

	
<b>Certificate of Compliance</b> RoHS Directive(2011/65/EU) Of the European Parliament and of the Council on the Restriction of the use of certain Hazardous Substances in Electrical and Electronic equipment (Pb, Cd, Hg, CrVI, PBBs, PBDEs, Phthalates) Certificate Number: CER21007003E      Issued Date: Aug 27, 2021	
Applicant:	Shenzhen Yuejiang Technology Co., Ltd.
Address:	Room 1003, Building 2, Chongwen Park, Nanshan iPark, No.3370, Liuxian Blvd, Fuguang Community, Taoyuan Street, Nanshan District, Shenzhen
Sample Name:	DOBOT MG400/Magician Pro
Sample Model:	DT-MG400-4R075-01, DT-MGPro-4R075-01
Sample Brand:	 DOBOT
This certificate of compliance is based on a test procedure or an evaluation of the above-mentioned product. This is to certify that the above-mentioned product is in compliance with the RoHS Directive (2011/65/EU) of the European Parliament. It is only valid in connection with the test report (Report Number: HLF21007003E)	
<b>RoHS</b>	Signed for and on Behalf of FTT  Technology Manager
Shenzhen FLION Testing Technologies Co., Ltd Gangzi industrial park, Furong industrial zone, Shajing, Bao'an, Shenzhen, China Tel:0755-27248885 Fax:0755-27480090 Website:Http://www.cnftt.com	

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<b>GACA</b>	Global Access Certification Australia Pty. Ltd. 960 Waverley Road, Wheelers Hill, Victoria 3150, Australia
<b>Proof of RCM Compliance</b>	
<b>Certificate No.: RCMP21632 001</b>	
Certificate holder:	Shenzhen Yuejiang Technology Co., Ltd. Room 1003, Building 2, Chongwen Park, Nanshan iPark, No. 3370, Liuxian Blvd, Fuguang Community, Taoyuan Street, Nanshan District, Shenzhen, Guangdong, China
Product Description:	DOBOT MG400 / DOBOT Magician Pro
Brand Name:	DOBOT
Model Number:	DT-MG400-4R075-01, DT-MGPro-4R075-01
Electrical Ratings:	48Vdc, 3.2A, 150W
<b><i>GACA approves that this Electrical Article as described above meets the requirements of the standard/regulation(s) listed below:</i></b>	
Complies with	Standard/regulation(s)
ACMA mandated	EN IEC 61000-6-4: 2019, EN IEC 61000-6-2: 2019, IEC 61000-3-2: 2014, IEC 61000-3-3: 2013+A1: 2017 (Test report: SZEM210800726801, SZEM210100007201)
Electrical Safety	IEC 60204-1: 2016, EN ISO 10218-1: 2011 (Test report: SZFS210800002501, SZFS210810002501)
EESS-National Equipment Database Registration	Above models is registered as level 1 product on <a href="https://equipment.era.gov.au/Registration/EquipmentSearch.aspx?atn=public">https://equipment.era.gov.au/Registration/EquipmentSearch.aspx?atn=public</a>
Date of Registered:	18 October 2021
Date of Expiry:	17 October 2026
 On behalf of GACA	

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